

# Product data sheet

✉ [sales@brannan.co.uk](mailto:sales@brannan.co.uk)

☎ +44 (0)1946 816600

🌐 [brannan.co.uk](http://brannan.co.uk)

## INFRARED EAR THERMOMETER

A quick and easy infrared thermometer suitable for all ages.  
Electronic thermometer with digital display and memory.  
Complete with protective bag and replacement probe covers.  
Especially suited with use with babies, young children and the elderly.

### Product features:

- Easy and safe
- Selectable °C or °F
- Supplied with battery, full instructions and protective bag

### Product specifications:

<b>Measuring temperature range:</b>	+35.5 to +42°C&F
<b>Displayed temperature range:</b>	+34 to +44°C&F
<b>Operating temperature range:</b>	+16 to +35°C (+60.8 to +95.0°F)
<b>Operating relative humidity:</b>	<80%RH
<b>Resolution:</b>	0.1°C&F
<b>Accuracy:</b>	+/-0.2°C between 35 and 42°C (+/-0.4°F between 95.9 and 107.6°F)
<b>Battery:</b>	2 x AAA
<b>CE, ROHS/WEEE:</b>	Approved
<b>Hazard information (SDS):</b>	See <a href="http://brannan.co.uk">brannan.co.uk</a> for information
<b>LCD dimensions:</b>	10mm x 20mm
<b>Product dimensions:</b>	Net: 140mm x 70mm x 40mm Gross: 173mm x 116mm x 49mm
<b>Product weight:</b>	Net: 67g (including batteries) Probe covers: 18g Display box weight: 75g Gross: 175g



Description	Barcode	Product No
Infrared ear thermometer	5 011405 114204	11/420/3



View me on  
the website

