

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk



SLOT FRONT LABORATORY THERMOMETER CASES

To protect the standard 6.5mm diameter laboratory thermometer from breakage, a variety of cases are available. These have a slot front for easy reading, a suspension ring at the top and holes in the bulb end to ensure quick response to temperature.

Product features:

- Suitable for 155mm, 305mm 355mm and 405mm thermometer lengths
- Available in brass or stainless steel
- Suspension ring at top of case
- Holes at bulb end of case for quick response
- Stainless steel scoop case also available for temperature measurement by sampling

Description	Product weight	Product no
155mm Brass slot front case	33g	47/732/0
305mm Brass slot front case	55g	47/432/0
355mm Brass slot front case	63g	47/832/0
405mm Brass slot front case	67g	47/932/0
305mm Stainless slot front case	61g	47/434/0
405mm Stainless slot front case	75g	47/934/0
305mm Stainless scoop slot front case	91.5g	47/435/0



View me on
the website

