

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

SPH012 HYDROMETER - SPECIAL PURPOSE HYDROMETER - URINOMETER

These hydrometers are 160mm long and measure the specific gravity of urine at 15.6°C (60°F). All hydrometers can be supplied with UKAS certification, certificate traceable to Recognised National Standard or Statement of Conformity, if required.

Product features:

- Measures the specific gravity of urine
- Adjusted at 15.6°C (60°F)
- 160mm long
- Made in UK
- Supplied in protective, resealable plastic tube

Product specifications:

Material: Soda glass, steel shot ballast, synthetic wax resin
Range: 1.000 to 1.060°SG
Accuracy: +/- 0.002°SG
Divisions: 0.001°SG
Bulb diameter: 16mm
Overall length: 160mm
Weight: 21g (including protective plastic tube)

Standard product numbers:

Range	Divisions	Product no
1.000 to 1.060°SG	0.001°SG	10/650/0



View me on
the website

