

# Brannan Catalogue

## Medical & Healthcare Range






- Personal Healthcare
- Environmental Monitoring
- Veterinary Instruments

**Measurement for Life.**



# S.Brannan & Sons



S. Brannan & Sons has manufactured and supplied thermometers, gauges and instrumentation since 1913. Over a century of tradition, combined with modern manufacturing methods at our facilities in the UK ensure unequalled service.

The Company has operated a quality management system since 1991 and currently operates within the scope of ISO 9001: 2015.

Complimenting our production facilities, the Company provides a calibration service operating its own UKAS accredited laboratory where, in addition to having new thermometers tested, customers may also have existing instruments re-calibrated.

Following the 2014 EU ban on mercury, a wide choice of electronic instruments and organic, non-hazardous filling materials, including our Patented LO-**tox**™ fillings, are available.

Products in this catalogue are generally supplied on Brannan standard packaging, or customer's own packaging (quantity dependant) and many are also available bulk packed.

Customers' requirements can be discussed with our Sales Team.

Products may also be viewed our website: [www.brannan.co.uk](http://www.brannan.co.uk)

Please note:

- All electronic instruments are supplied with batteries and full instructions.
- All dimensions are for unpacked products
- All information is correct at time of printing but our policy is one of continuous product development. We reserve the right to change specification without notice.



## EMAIL

[sales@brannan.co.uk](mailto:sales@brannan.co.uk)

## PHONE

+ 44 (0)1946 816600

## WEBSITE

[www.brannan.co.uk](http://www.brannan.co.uk)

## FAX

+ 44 (0)1946 816625

## ADDRESS

S.Brannan & Sons  
Leconfield Industrial Estate  
Cleator Moor  
Cumbria  
UK  
CA25 5QE

[www.brannan.co.uk](http://www.brannan.co.uk)

**Measurement for Life.**



**Brannan** offer an extensive range of healthcare thermometers for professional or personal use. The following products are our easy to use thermometers suitable for use in the home to give quick and clear readings. All our digital instruments are CE, WEEE and RoHS compliant and are supplied complete with instructions and batteries (where appropriate).



## 1. Standard Oral Clinical Thermometer

Rigid probe thermometer supplied in protective case.

Suitable for use with disposable probe covers

(Ref: 11/066/0).

Range: +32 to +42.9°C

Dimensions: 133 x 24 x 14mm

LCD Dimensions: 6 x 18mm

11/064/2

**BEST  
SELLER**

UKAS  
Calibration  
Available

## 2. Flexi-head Oral Clinical Thermometer with Large Display

Flexible probe clinical thermometer for gentle temperature measurement. Suitable for use with disposable probe covers

(Ref: 11/066/0).

Range: +32 to +41.9°C

Dimensions: 130 x 22 x 12mm

LCD Dimensions: 10 x 24mm

11/067/2



**NEW**

## 3. Waterproof, Flexi-head Oral Clinical Thermometer with Disposable Probe Covers

Advanced, waterproof, flexible probe clinical thermometer with high speed response reading and protective case. Supplied with 25pcs disposable probe covers for hygienic use of thermometer.

Range: +32 to +42.9°C


Dimensions: 133 x 22 x 12mm

LCD Dimensions: 7 x 20mm

11/065/2



Description	Range	Product Dimensions	Product Number
1. Standard Clinical Thermometer	+32 to +42.9°C	133 x 24 x 14mm	11/064/2
2. Large LCD Flexi-Head Clinical Thermometer	+32 to +41.9°C	130 x 24 x 14mm	11/067/2
3. Waterproof Flexi-Head Clinical Thermometer	+32 to +42.9°C	133 x 22 x 12mm	11/065/2
Pack 25 Disposable Probe Covers	n/a	130 x 50 x 30mm	11/066/0



## Ovulation Kit

By taking your temperature in a controlled and consistent manner each morning and plotting the results on a chart, it is possible to follow the progress of your menstrual cycle.

Oral clinical thermometer supplied with book of recording charts and instructions for daily readings.

Dimensions: 210 x 150 x 15mm

11/020/3

NB - This is not a contraceptive device.

To interpret the results you must seek the advice of your doctor.





**Market leading  
Used by UN organisations**

**Advanced Body & Surface  
Temperature Non-Contact,  
Infrared Thermometer**

Non-contact and hygienic infrared thermometer designed to take the body temperature of a person regardless of room temperature. Ideally suited for use in schools, businesses or other environments where a large number of people are coming into contact with one another. Features include selectable body or surface temperature and switchable °C/°F readout, alarm value setting, trigger lock and data hold.

Memory Recall: 32 measurements  
Body Temp Range: +32 to +42.5°C&F  
Surface Temp Range: 0 to +60°C&F  
Dimensions: 149 x 77 x 43mm  
LCD Dimensions: 28 x 24mm  
Pip Code: 347-3733  
**11/400/3**

Statement of  
Conformity  
Available

Traceable  
Calibration  
Available



**Infrared Ear Thermometer**

Easy to use electronic thermometer with digital display and memory of last 10 readings. Complete with protective bag and replacement probe covers for hygienic use of ear thermometer.


Range: +35.5 to +42°C&F  
Dimensions: 140 x 70 x 90mm  
LCD Dimensions: 10 x 20mm  
Pip Code: 289-2305  
**11/420/3**



**Infrared Ear Thermometer  
Replacement Probe Covers (20pcs)**

Boxed Dimensions: 76 x 32 x 30mm  
Pip Code: 329-4303  
**11/421/0**

**BEST  
SELLER**



**Forehead Strip Temperature Indicator**

This quick and safe forehead strip thermometer uses heat sensitive liquid encased in plastic film. Ideal for taking the temperature of small children. Comes complete with protective case and full instructions.

Ideal for use in First Aid Kit.  
Range: +35 to +40°C&F  
Dimensions: 39 x 99 x 10mm  
PIP Code: 047-6143  
**11/465/3**



## Automatic, Non Contact, Rapid Test, Infrared Temperature Scanner

Introducing the automatic temperature scanner specifically designed to assist in current battle against COVID-19 and other infections.

The robust thermometer can be wall mounted or positioned on a tripod, situated in a position where a high number of people are entering a premises and capable of taking multiple temperature readings within a short period of time. Ideally suited for the entrance to schools, office buildings, retail shops, gyms, etc, and eliminates the need for a member of staff to approach visitors.

**Please note:** This is not a medical device and should not be used for diagnostic purposes. We recommend the use of a medical body thermometer to confirm high temperature readings.



### Product Features:

- Automatic temperature measurement - hygienic, non contact
- Distinctive audible sound and red flashing light alert when high temp detected
- Connects to PC for remote reading and data storage
- USB cable power supply (suitable for use with rechargeable lithium battery\*)
- Quick response reading (0.5s)
- Long measuring distance
- High accuracy – 0.2°C
- Wall fixture or optional tripod stand\*

**Temperature Scanner : 14/100/0**


**Scanner Tripod : 14/100/1**

**Rechargeable Battery: 14/100/2**

\*Additional charge applies

## Hypothermia Room Thermometers

The elderly are particularly susceptible to fluctuations in temperature. It is therefore important to maintain a balanced environment for them. These thermometers can be supplied with a Certificate of Conformity or a Calibration Certificate traceable to National Standards.



### Hypothermia Wall Thermometers

Choice of 240mm, 215mm or 145mm wall thermometers with zones to help determine the most comfortable and economical temperature settings for your home.

240 x 65 x 8mm - **14/867/2**

145 x 34 x 7mm - **14/862/2**

Range: -20 to +50°C only

215 x 40 x 8mm - **14/866/2**

Statement of Conformity Available


### 65mm Dial Thermometer on Pyramid card

Range: +20 to +120°F&C

Dimensions: 150 x 150 x 80mm

Pip Code: 320-1829

**14/865/3**



### Hanging Dial Hypothermia Thermometer

65mm white dial with hook for hanging off shelf.

Large zones for easy reading.

Range: -10 to +50°C only

Dimensions: 85 x 65 x 12mm

Pip Code: 401-9311

**14/868/2**



### Medicine Storage Monitoring



#### Vaccine Thermometer

Vertical, hanging vaccine thermometer with zones to indicate the safe storage temperature for vaccines and medicines.

Range: -30 to +50°C only

Dimensions: 144 x 18 x 6mm

Pip Code: 320-1779

22/992/2

Statement of  
Conformity  
Available

Traceable  
Calibration  
Available

#### Blood Bank Bottle Thermometer

This unique combination of short scale thermometer fitted in a sealed bottle will smooth out temperature fluctuations caused by the operation of thermostats or opening the storage cabinet.

Range: -5 to +20°C only

Dimensions: 145 x 35mm dia

Pip Code: 320-1787

22/100/2



### Bath Water Monitoring

Brannan offer a range of bath thermometers to assist in obtaining the correct temperature of bath water, therefore avoiding the risk of scalding or chilling. This is invaluable for babies or elderly people.

#### Scoop Bath Thermometer

An alternative method of checking bath water temperature by sampling. Dual scale with zones indicating recommended temperatures.

Range: -25 to +50°C&F

Dimensions: 230 x 45 x 35mm

Pip Code: 031-3601

18/126/3



#### Horizontal Bath Thermometers

140mm white bath thermometers with zones showing recommended temperatures. This thermometer is ideally suited to customer own branding.

Available in standard black, dual scale design or fun and colourful Celsius only design.

Range: -10 to +50°C&F

Dimensions: 140 x 25 x 10mm

**Standard Design: 18/132/3 (°C&F)**

Pip Code: 031-3619

**Coloured Design: 18/140/2 (°C only)**



#### Baby Bath & Food Thermometer

Pocket test thermometers supplied in protective cases for checking temperature of baby's bath or food. Thermometer marked with zones to indicate safe temperature of water.

Range: 0 to +50°C&F

Dimensions: 160 x 20 x 16mm

Pip Code: 032-6942

18/348/3



Traceable  
Calibration  
Available

Statement of  
Conformity  
Available

#### Fish Shaped Bath Thermometer

A quick and fun way to check bath water temperature to avoid the risk of scalding or chilling. Showing recommended bath temperatures.

Range: -6 to +40°C only

Dimensions: 148 x 56 x 14mm

Pip Code: 291-2624

18/420/2





## Brooder Thermometer

Correct temperature control is essential for successful hatching. These thermometers will assist with the adjustment of thermostats and ventilation, helping to achieve optimum conditions during incubation. Choice of plastic or wooden brooder thermometers with easy to read scale and zones marked in weeks to show the various temperature levels required for brooding chicks.



Traceable  
Calibration  
Available

Statement of  
Conformity  
Available

### Plastic Brooder Thermometer

150mm short wall thermometer with week zones to indicate temperature levels required for brooding chicks. This thermometer is ideally suited for customers own branding or design.

Range: -10 to +50°C only  
Dimensions: 145 x 34 x 7mm  
20/155/2

## Digital Veterinary Thermometers

Easy to use electronic thermometers with chart showing Hobday's specification (the normal body temperature for a variety of animals). Standard, rigid probe thermometer or advanced, waterproof, flexi-head probe thermometer available, both supplied with protective case.

### Standard Vets Clinical Thermometer

Rigid probe thermometer supplied in protective case.

Range: +32 to +42.9°C

Dimensions: 133 x 24 x 14mm

LCD Dimensions: 6 x 18mm

11/015/2

**BEST  
SELLER**



### Digital Veterinary Thermometer

Hobday's Specification	Man	Cat	Horse	Dog	Cow	Pig	Sheep	Fowl
°C	37	38	38	38.8	38.9	39.5	40	41.7

### Waterproof, Flexi-head Oral Clinical Thermometer

Advanced, waterproof, flexible probe clinical thermometer with high speed response reading and protective case.

Range: +32 to +42.9°C

Dimensions: 133 x 22 x 12mm

LCD Dimensions: 7 x 20mm

11/035/2

**Terrarium Instruments**


Range of instruments which can be fixed inside a climate controlled chamber such as a terrarium for environmental monitoring.

**Mini Dial Thermometer & Hygrometers**


White dial instruments with green scales and red pointers. Thermometer, hygrometer and twin thermometer and hygrometer available with self adhesive pad for fixing inside terrarium. Particularly suited to customer own branding.

Ranges: -10 to +50°C only or °C&F  
0 to 100%RH

Mini Dial Description	Range	Product Dimensions	Product Number
Thermometer	-10 to +50°C only	65 x 65 x 10mm	30/440/2
	-10 to +50°C&F	65 x 65 x 10mm	30/440/3
Hygrometer	0 to 100%RH	65 x 65 x 10mm	30/441/0
Twin Thermometer & Hygrometer	°C only and %RH	130 x 65 x 10mm	30/442/2
	°C&F and %RH	130 x 65 x 10mm	30/442/3

**Pool & Pond Thermometers**

The monitoring of environmental conditions is important in keeping specialised or exotic fish. These thermometers will monitor the temperature of pond water for carp or other fish, or the temperature of your pool.

**Floating Pool Thermometer**

Blue float supports a white thermometer with blue temperature scale.  
Range: -10 to +50°C&F  
Dimensions: 175 x 60 x 70mm  
**18/416/3**

**Scoop Pool & Pond Thermometer**  
230mm blue scoop pool and pond thermometer with white temperature scale. An alternative method of checking water temperature by sampling.

Range: -10 to +50°C&F  
Dimensions: 245 x 46 x 30mm  
**18/410/3**

Statement of  
Conformity  
Available

Traceable  
Calibration  
Available





S. Brannan & Sons Ltd.  
Cleator Moor, Cumbria CA25 5QE, UK

Tel: +44 (0)1946 816600

Fax: +44 (0)1946 816625

Email: [sales@brannan.co.uk](mailto:sales@brannan.co.uk)



[facebook.com/brannanthermometersandinstrumentation](https://facebook.com/brannanthermometersandinstrumentation)



[plus.google.com/114245952167321950392/posts](https://plus.google.com/114245952167321950392/posts)



[twitter.com/Brannan\\_UK](https://twitter.com/Brannan_UK)

S/290018.Rev1.2021

[www.brannan.co.uk](http://www.brannan.co.uk)

**Measurement for Life.**

