

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

OVEN BATCH SAMPLING THERMOMETER

These precision test thermometers are sealed into a 30ml polypropylene bottle containing a heat transfer medium and are ideally suited to monitor the temperature of samples in controlled storage. The 'bottle' base of the assembly maintains a core temperature similar to that of the samples being monitored and is unaffected by short term fluctuations caused by thermostatic variations or through the opening of cabinet doors.

Brannan's unique **LO-tox™** blue filling eliminates the risk of mercury contamination. Supplied with Batch Sampling Thermometer Statement of Accuracy.

Product features:

- Non hazardous **LO-tox™** filling
- Temperature measurement unaffected by short term fluctuations
- Supplied with statement of accuracy
- Made in England

Product specifications:

Temperature range:	+20 to +130°C
Accuracy:	+/-2.0°C
Divisions:	1.0°C
CE/ROHS/WEEE:	N/A
Hazard Information (SDS):	See brannan.co.uk for information
Product Material:	Bottle: Polypropylene Bottle filling: Aluminium oxide crystals Glass tube filling: LO-tox™ filling
Product Weight:	Net: 78g Gross: 91g
Product Dimensions :	Net: ø 35mm x 145mm Display Box: 13g Gross: 153mm x 37mm x 37mm



Description	Barcode	Product No
Oven Batch Sampling Thermometer	5 011405 171023	17/102/2

View me on
the website

