

Product data sheet

✉ sales@brannan.co.uk

☎ +44 (0)1946 816600

🌐 brannan.co.uk

FOREHEAD STRIP TEMPERATURE INDICATOR

Quick and safe, forehead strip thermometer using heat sensitive liquid encased in plastic film. Ideal for small children or for adding into first aid kit.

Product features:

- Simple to use and easy to read
- Quick and safe
- Supplied with storage sleeve
- Made in EU

Product specifications:

Temperature range:	+37 to +40°C&F
Indicators:	Every 1.0°C/°F
Accuracy:	+/- 1.0°C
Storage & operating conditions:	16 to 40°C (60 to 104°F) at 15 to 95% RH*
CE & RoHS/WEEE:	Compliant
Hazard Information (SDS):	See brannan.co.uk for information
Product dimensions:	Strip: 89mm x 19mm Sleeve: 111mm x 79mm x 3mm
Product weight:	Strip: 0.5g Gross: 9g


*Note:

After use, store thermometer in the protective sleeve and keep away from Ultraviolet/sunlight, dampness, or heat.

The thermometer can be cleaned by gently wiping with a soft cloth or tissue dipped in water.

Do not immerse thermometer in liquids

Optimal storage and usage between 16 to 27°C (60 to 80°F) at 40 to 60% RH.



Description	Barcode	Product No
Forehead temperature indicator	5 011405 114662	11/466/9

